

WellisAir Disinfection Purifier for Hospital Safety

Wellis removes microbial contaminants both airborne and on the surface because WellisAir emits hydroxyl radicals (OH) in all directions.



5-minute installation and easy maintenance

- Wide coverage (100m³ - 300m³)
- Air quality Display LED
- Mode control (day/night)
- Cartridge sensing and display

Simple controls and low-cost maintenance



Automatic operation

Automatically purifies airborne and surface contaminants



Night Mode

Minimal noise, quiet functioning. Low brightness for comfortable sleep



Cartridge Alarm

Lets you know in advance time to replace cartridge



Alert

Notifies in case of empty cartridge or when improper inserting causes a stop in function

Laboratory test results

Items	Result		
	Initial density	4 hours later (CFU/mL)	Bacteria reduction (%)
Antibacterial test for E. coli	1.7x10 ⁴	<10	99.9
Antibacterial test for Pseudomonas aeruginosa	1.9x10 ⁴	<10	99.9
Antibacterial test for Staphylococcus aureus	1.4x10 ⁴	<10	99.9
Antibacterial test for Salmonella	1.6x10 ⁴	<10	99.9
Antibacterial test for Pneumonia	2.0x10 ⁴	<10	99.9
Antibacterial test for MRSA bacteria	1.2x10 ⁴	<10	99.9

Lab results shows a 99.9% reduction in bacteria pathogens during an antibacterial test run to verify the existence of bacterial contaminants on surfaces such as clothing and furniture.

9LRLSSPPRRL
 I! IRLSSPNRRL
 IRLRLSRRL
 IPS! PNRRL
 7000fi 6635839

L



WellisAir